


# Search for a Freezer

Pg. 1 of 1


An EIDSS user with the appropriate permissions can locate freezers in the database. This quick guide presents the steps to search for a freezer.

1. From the **EIDSS Dashboard**, choose **FREEZERS** from the **LABORATORY** drop-down menu. The **Freezer List** screen appears.



2. Click **Advanced Search**. The **Advanced Search** window appears.


**Tip 1:** To perform a basic search, enter a full or partial **Freezer Name** in the text box and click the **Search** icon.



3. Enter as much information as possible about the freezer for which you are searching.

**Tip 2:** To delete entered information and close the window, click **Cancel**.

**Tip 3:** To delete entered information, click **Clear**.



4. Click **Search**. The **Freezer List** screen refreshes displaying the freezers that match your search criteria.