


Export a Layout


Pg. 1 of 2

An EIDSS user with the appropriate permissions can export layouts as spreadsheets or other document file formats from the Analysis, Visualization, and Reporting (AVR) application. This quick guide presents the steps to export a layout.

1. From the AVR **Query and Layout Tree**, highlight the layout you wish to export.



2. Choose **Report** from the **Export** drop-down menu. The format options appear.
3. Choose the desired file format. A **Save As** window appears.



Export a Layout

Pg. 2 of 2

4. Navigate to the location where you wish to save the layout.
5. Enter the **File name**. **Note:** The **Save as type** field is auto-populated based on the chosen format.
6. Click **Save**. The layout is exported to the chosen location.

