


Accession a Sample


Pg. 1 of 2

An EIDSS user with the appropriate permissions can add samples to an existing collection in the database. This quick guide presents the steps to accession a sample.

 From the EIDSS Dashboard, choose SAMPLES from the LABORATORY drop-down menu. The Laboratory screen appears with the SAMPLES tab open displaying all samples either unaccessioned (highlighted in yellow) or accessioned and not assigned a test.



Check the box of each sample you wish to accession. The Laboratory screen refreshes with the Accession button enabled.





Accession a Sample

Pg. 2 of 2

- 3. Click **Accession**. A pop-up appears offering additional fields.
- Choose a Sample Status from the drop-down menu. Note: If multiple samples were selected, the status chosen will apply to all.
 Note: If Accept in Poor Condition or Reject is chosen, a red asterisk will appear next to the corresponding Comment icon in the grid indicating a mandatory field.
- Check the Print Barcode? box as needed. The Laboratory screen refreshes with the modified sample bolded and the Save button enabled. Note: If the box is checked, click Print on the Bar Code Labels window that appears after saving.
- Click Save. A pop-up appears confirming successful accessioning of the sample.

