


Assign a Test


Pg. 1 of 2

An EIDSS user with the appropriate permissions can assign tests to accessioned samples within the TESTING tab in the database. This quick guide presents the steps to assign a test.

1. From the **EIDSS Dashboard**, choose **TESTING** from the **LABORATORY** drop-down menu. The **Laboratory** screen appears with the **TESTING** tab open.



2. Enter a search parameter for the sample to which you wish to assign a test. **Note:** Searching may not be necessary if the desired sample appears in the grid.



3. Click the **Search** icon. The **Laboratory** screen refreshes displaying the **Search Results**.

Tip 1: If the desired sample does not appear, click **Advanced Search**, enter as much information as possible, and click **Search**.



4. Check the box of the desired sample.
5. Click **Assign Test**. The **Assign Test** window appears.

Assign a Test

Pg. 2 of 2

6. Choose a **Test Name** from the drop-down menu. **Note:** Uncheck the **Filter Test Name by Disease** box if you do not wish to filter the drop-down menus by the test names and categories associated with the disease.
7. Click the **Add** icon. **Note:** Repeat steps 6 and 7 to assign additional tests to the sample.

Tip 2: To delete entered information and close the window, click **Cancel**.

8. Click **Save**. The **Laboratory** screen refreshes with the modified sample bolded and the **Save** button enabled.
9. Click **Save**. A pop-up appears confirming successful assignment of the test.

	Report/Session ID	Patient/Farm Owner	Sample Type	Disease	Test Name
<input type="checkbox"/>	HUM0000190081	Smith Rodney	Aspirate - Lymph	Acute Respiratory Infection	Dengue test
<input type="checkbox"/>	HUM0000190080	Smith Rodney	Aspirate - Nasal	Acute Respiratory Infection	Dengue test
<input type="checkbox"/>	HUMTBTBZ190497	Williams Jeff	Blood - anticoagulated whole blood	Brucellosis - Suis	Culture