

Configure the Veterinary Prophylactic Measure Matrix


Pg. 1 of 4

An EIDSS user with the appropriate permissions can manage the list of prophylactic measures in the database. This quick guide presents the steps to complete the following tasks:

- Create a New Veterinary Prophylactic Measure Matrix Version
- Edit an Existing Veterinary Prophylactic Measure Matrix Version

Create a New Matrix Version

1. From the **EIDSS Dashboard**, choose **VETERINARY PROPHYLACTIC MEASURE MATRIX** from the **CONFIGURATIONS** drop-down menu. The **Vet Prophylactic Measure Matrix** screen appears.



2. If creating a new matrix from a previous version, choose the desired version from the **Select Version** drop-down menu. **Note:** You can choose either an active or inactive version; only one version can be active at a time. **Note:** If no previous matrix versions exist, proceed to step 6.

The screenshot shows the 'Matrix Version' form. The 'Select Version' dropdown menu is open, showing 'LJM Copy' as the selected version. The form also includes fields for 'Version Name' and 'Activation Date'. Below the dropdown, there are buttons for 'Activate Matrix Version', 'New Matrix Version', and 'Delete Matrix Version'. The 'Vet Prophylactic Measure Matrix' section is visible below the form.

3. Click **New Matrix Version**. A confirmation dialog box appears.

The screenshot shows a 'Confirmation' dialog box with the text 'Do you want to save changes?'. There are 'Yes' and 'No' buttons. The dialog box is overlaid on the 'Matrix Version' form.

4. Click **Yes**. A pop-up appears.

The screenshot shows a confirmation dialog box with the text 'The record is saved successfully.' and an 'OK' button. The dialog box is overlaid on the 'Matrix Version' form.

5. Click **OK**. The screen refreshes displaying a copy of the desired version.

Configure the Veterinary Prophylactic Measure Matrix

Pg. 2 of 4

6. Enter the **Version Name**.

7. Enter an **Activation Date**.

Tip 1: To navigate to the desired date, click the **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

8. Click **Add**. A new row appears.

9. Choose a **Measure Type** from the drop-down menu. The **Measure Code** auto-populates.

10. Choose a **Species** from the drop-down menu.

11. Choose a **Disease** from the drop-down menu. The **OIE Code** auto-populates..
Note: Repeat steps 8–11 to add additional rows.

Tip 2: To add a new option to a drop-down menu, click the appropriate **Add** icon, complete the window that appears, and click **Save**.

Tip 3: To change the display order of rows in the grid, highlight the desired row and drag it to the appropriate location.

Tip 4: To delete entered information, click **Cancel**.

12. Click **Save**. A confirmation dialog box appears.

The screenshot shows the 'Matrix Version' form. At the top, there are fields for 'Version Name' (Prophylactic_Training), 'Activation Date' (6/5/2019), and 'Select Version' (LJM Copy Copy). Below these are buttons for 'Activate Matrix Version', 'New Matrix Version', and 'Delete Matrix Version'. The main section is titled 'Vet Prophylactic Measure Matrix'. It features a table with columns: Order, Measure Type, Measure Code, Species, Diagnosis, and OE Code. The first row (Order 1) contains: ABCGGG1234, ABCGGG1234, Cattle, Avian chlamydiosis. Below the table, there are navigation buttons: 'Add', 'Delete Matrix Record', 'Save', 'OK', and 'Cancel'. A callout 'T1' points to the 'Add' button.

This screenshot shows the same 'Matrix Version' form, but now with two rows in the table. The second row (Order 2) contains: GGGG1, GGGG1, Dog, Brucellosis - Unspecified. Callouts are present: 'T2' points to the 'Measure Type' dropdown of the second row; '9' points to the 'Add' icon next to the second row; '10' points to the 'Species' dropdown of the second row; '11' points to the 'Diagnosis' dropdown of the second row; '12' points to the 'Save' button; and 'T4' points to the 'Cancel' button. The table footer shows 'Showing 1 to 2 of 2 entries'.


Configure the Veterinary Prophylactic Measure Matrix

Pg. 3 of 4

13. Click **Yes**. A pop-up appears confirming successful saving of the record.



14. Click **OK**.



15. Click **Activate Matrix Version**. A pop-up appears confirming successful configuration of the matrix.

A screenshot of the 'Matrix Version' configuration form. The form has a title bar 'Matrix Version' and contains the following fields: 'Version Name' (text input with 'Prophylactic_Training'), 'Activation Date' (date picker with '6/5/2019'), and 'Select Version' (dropdown menu with 'LJM Copy Copy'). Below these fields are three buttons: 'Activate Matrix Version', 'New Matrix Version', and 'Delete Matrix Version'. An orange circle with the number '15' points to the 'Activate Matrix Version' button. Below the buttons is a section titled 'Vet Prophylactic Measure Matrix'. This section includes a 'Show' dropdown set to '10' and 'entries'. Below this is a table with columns: 'Order', 'Measure Type', 'Measure Code', 'Species', 'Diagnosis', and 'OE Code'. The table contains one entry with 'Order' 1, 'Measure Type' 'ABCGGG1234', 'Measure Code' 'ABCGGG1234', 'Species' 'Cattle', and 'Diagnosis' 'Avian chlamydiosis'. Below the table is a 'Showing 1 to 1 of 1 entries' message and 'Previous', '1', 'Next' navigation buttons. At the bottom of the form are buttons for 'Add', 'Delete Matrix Record', 'Save', 'OK', and 'Cancel'.

Configure the Veterinary Prophylactic Measure Matrix

Pg. 4 of 4

Edit an Existing Matrix Version

- Choose the version you wish to edit from the **Select Version** drop-down menu. The version appears in the grid.
Note: In the **Select Version** drop-down menu, the activated version has a red checkmark, inactive versions have green checkmarks, and never-activated versions have no checkmarks.

Tip 1: To activate a never-activated version, choose it from the **Select Version** drop-down menu and click **Activate Matrix Version**.

Tip 2: To delete an inactive or never-activated version, choose it from the **Select Version** drop-down menu, click **Delete Matrix Version**, and click **Yes** in the dialog box that appears.

- Revise the **Version Name** and **Activation Date** as needed. **Note:** The grid is editable only for never-activated versions.

Tip 3: To add a row to a never-activated version, click **Add**, complete the fields that appear, and click **Save**.

Tip 4: To delete a row from a never-activated version, highlight it in the grid, click **Delete Matrix Record**, and click **Yes** in the dialog box that appears.

- Click **Save**. A confirmation dialog box appears.
- Click **Yes**. A pop-up appears confirming successful configuration of the matrix.

Matrix Version

Version Name: LJM Activation Date: 05/30/2019 Select Version: **LJM** (ACTIVATION DATE: 05/30/2019)

Buttons: Activate Matrix Version, New Matrix Version, Delete Matrix Version

Vet Prophylactic Measure Matrix

Show 10 entries

Order	Measure Type	Measure Code	Species	Diagnosis	OE Code
1	Vaccination	vac	Mule	Enzootic bovine leucosis	
2	GGGG1	GGGG1	Boar	African Swine Fever	A120
3	GGGG1222	GGGG1222	Goat	Avian chlamydiosis	
4	Test1234	Test1234	Horse	Brucellosis - Melitensis	
5	Test33	Test33	WBEAST2	Brucellosis - Unspecified	

Matrix Version

Version Name: LJM Copy Test Activation Date: 6/5/2019 Select Version: **LJM** (ACTIVATION DATE: 05/30/2019)

Buttons: Activate Matrix Version, New Matrix Version, Delete Matrix Version

Vet Prophylactic Measure Matrix

Show 10 entries

Order	Measure Type	Measure Code	Species	Diagnosis	OE Code
1	Vaccination	vac	Mule	Enzootic bovine leucosis	
2	GGGG1	GGGG1	Boar	African Swine Fever	A120
10	ABCGGG1234	ABCGGG1234	Aardvark	Bovine tuberculosis	

Showing 1 to 10 of 15 entries

Buttons: Add, Delete Matrix Record, **Save**, OK, Cancel

Confirmation

Do you want to save changes?

Buttons: Yes, No