


Edit a Transfer

An EIDSS user with the appropriate permissions can revise transfer details in the database. This quick guide presents the steps to edit a transfer.


 From the EIDSS Dashboard, choose TRANSFERRED from the LABORATORY drop-down menu. The Laboratory screen appears with the TRANSFERRED tab open.




- Enter a search parameter for the transferred sample you wish to edit.
 Note: Searching may not be necessary if the desired sample appears in the grid.
- Click the Search icon. The Laboratory screen refreshes displaying the Search Results.



 Click the Edit icon corresponding to the desired sample. The Sample/Test Details window appears displaying the Transfer Details section.







Edit a Transfer

5. Make any necessary revisions.

Tip 2: To navigate to the desired date, click the Calendar icon and then the Select Month >Select Year > Select Decade bars at the top of the calendar.

Tip 3: To delete entered information and close the window, click **Cancel**.

- Click Save. The Laboratory screen refreshes with the modified sample bolded and the Save button enabled.
- 7. Review the revised information for accuracy.
- 8. Click **Save.** A pop-up appears confirming successful modification of the transfer information.

