


Configure the Test–Test Results Matrix

Pg. 1 of 1

An EIDSS user with the appropriate permissions can manage the lists of test results as specified by test name in the database. This quick guide presents the steps to configure the Test–Test Results matrix.

- From the **EIDSS Dashboard**, choose **TEST – TEST RESULT MATRIX** from the **CONFIGURATIONS** drop-down menu. The **Test–Test Results Matrix** screen appears.



- Select the radio button corresponding to the desired test result type.
- Choose a **Test Name** from the drop-down menu. Any existing **Test Results** appear in the grid.

The screenshot shows the 'Test-Test Results Matrix' form. In the 'Test Results Relations' section, the 'Lab Test' radio button is selected. A dropdown menu for 'Test Name' is open, showing 'Acetic agar' with a red circle and a '2' callout. Below the dropdown is a table with columns 'Row', 'Test Result', and 'Indicative'. The table contains three rows: 'Growth neg', 'Growth pos', and 'Indeterminate'. A red circle and a '3' callout are on the 'Test Result' column. A red circle and a '4' callout are on the 'Add' button. A red circle and a 'T1' callout are on the 'Test Name' dropdown.

Tip 1: To add a new option, click the **Add** icon, complete the window that appears, and click **Submit**.

Tip 2: To delete a row, click the corresponding **Delete** icon and click **Yes** on the dialog box that appears.

- Click **Add**. The **Test–Test Results Matrix** window appears.
- Choose a **Test Result** from the drop-down menu.

Tip 3: To mark a **Test Result** as indicating detection of the disease, check the **Indicative** box; the **Indicative** box can be checked for multiple **Test Results**.

Tip 4: To delete entered information and close the window, click **Cancel**.

The screenshot shows the 'Test-Test Results Matrix' form. In the 'Test Result' section, the 'Positive' dropdown menu is open with a red circle and a '5' callout. Below the dropdown is a checkbox labeled 'Indicative' which is checked with a red circle and a 'T3' callout. At the bottom of the form are 'Submit' and 'Cancel' buttons. A red circle and a '6' callout are on the 'Submit' button. A red circle and a 'T4' callout are on the 'Cancel' button. A red circle and a 'T1' callout are on the 'Test Name' dropdown.

- Click **Submit**. A pop-up appears confirming successful configuration of the matrix.