


Delete a Test Record


Pg. 1 of 2

An EIDSS user with the appropriate permissions can remove records from the database. This quick guide presents the steps to delete a test record.

1. From the **EIDSS Dashboard**, choose **TESTING** from the **LABORATORY** drop-down menu. The **Laboratory** screen appears with the **TESTING** tab open.



2. Enter a search parameter for the sample with the test you wish to delete. **Note:** Searching may not be necessary if the desired sample appears in the grid.



3. Click the **Search** icon. The **Laboratory** screen refreshes displaying the **Search Results**.

Tip: If the desired sample does not appear, click **Advanced Search**, enter as much information as possible, and click **Search**.




4. Check the box of each desired sample. **Note:** Only accessioned samples with a **Test Name** and a **Test Status** of **Preliminary** or **In Progress** can be deleted.
5. Click the **Samples Menu** icon. The **Samples Menu** appears.

Delete a Test Record

Pg. 2 of 2

- Click **Lab Record Deletion > Delete Test**. The **Laboratory** screen refreshes with the modified sample bolded and the **Save** button enabled.



- Click **Save**. The test is moved to the **APPROVALS** tab.

