


# Configure the Disease–Sample Type Matrix

Pg. 1 of 1


An EIDSS user with the appropriate permissions can manage the list of sample types specified by disease in the database. This quick guide presents the steps to configure the Disease–Sample Type matrix.

- From the **EIDSS Dashboard**, choose **DISEASE – SAMPLE TYPE MATRIX** from the **CONFIGURATIONS** drop-down menu. The **Disease–Sample Type Matrix** screen appears displaying any existing **Sample Types**.



**Tip 1:** To delete a row, click the corresponding **Delete** icon and click **Yes** on the dialog box that appears.


- Click **Add**. The **Disease–Sample Type Matrix** window appears.



- Choose a **Disease** from the drop-down menu.
- Choose a **Sample Type** from the drop-down menu.

**Tip 2:** To add a new option, click the appropriate **Add** icon, complete the window that appears, and click **Submit**.

**Tip 3:** To delete entered information and close the window, click **Cancel**.



- Click **Submit**. A pop-up appears confirming successful configuration of the matrix.