


Configure the Vector Type – Sample Type Matrix

Pg. 1 of 1

An EIDSS user with the appropriate permissions can manage lists of sample types specified by vector type in the database. This quick guide presents the steps to configure the Vector Type – Sample Type matrix.


1. From the **EIDSS Dashboard**, choose **VECTOR TYPE – SAMPLE MATRIX** from the **CONFIGURATIONS** drop-down menu. The **Vector Type – Sample Type Matrix** screen appears.



2. Choose a **Vector Type** from the drop-down menu. Any existing **Sample Types** appear in the grid.


Tip 1: To add a new option, click the **Add** icon, complete the window that appears, and click **Submit**.

Tip 2: To delete a row, click the corresponding **Delete** icon and click **Yes** on the dialog box that appears.



3. Click **Add**. The **Vector Type – Sample Type Matrix** window appears.
4. Choose a **Sample Type** from the drop-down menu.

Tip 3: To delete entered information and close the window, click **Cancel**.



5. Click **Submit**. A pop-up appears confirming successful configuration of the matrix.