


# Edit a Test


Pg. 1 of 2

An EIDSS user with the appropriate permissions can revise test details in the database. This quick guide presents the steps to edit a test.

- From the **EIDSS Dashboard**, choose **TESTING** from the **LABORATORY** drop-down menu. The **Laboratory** screen appears with the **TESTING** tab open.



- Enter a search parameter for the sample with the test you wish to edit.  
**Note:** Searching may not be necessary if the desired sample appears in the grid.



- Click the **Search** icon. The **Laboratory** screen refreshes displaying the **Search Results**.

**Tip 1:** If the desired sample does not appear, click **Advanced Search**, enter as much information as possible, and click **Search**.



- Click the **Edit** icon corresponding to the desired sample. The **Sample/Test Details** window appears displaying the **Test Details** section.

## Edit a Test

Pg. 2 of 2

5. Make any necessary revisions.
6. Use the **Show/Hide** icons to access and revise other sections as needed.

**Tip 2:** To delete entered information and close the window, click **Cancel**.

7. Click **Save**. The **Laboratory** screen refreshes with the modified sample bolded and the **Save** button enabled.

8. Review the revised information for accuracy.
9. Click **Save**. A pop-up appears confirming successful modification of the test.

Report/Session ID	Patient/Farm Owner	Sample Type	Disease	Test
<b>HUM0000190081</b>	<b>Smith Rodney</b>	<b>Aspirate - Lymph</b>	<b>Acute Respiratory Infection</b>	<b>Dengue test</b>
HUM0000190080	Smith Rodney	Aspirate - Nasal	Acute Respiratory Infection	Dengue test
HUMTBTBZ190497	Williams Jeff	Blood - anticoagulated whole blood	Brucellosis - Suis	Culture