


# Delete a Vector Surveillance Session

Pg. 1 of 2

An EIDSS user with the appropriate permissions can remove surveillance sessions from the database. This quick guide presents the steps to delete a vector surveillance session.

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears. **Note:** All field tests and samples within detailed collections; detailed collections; aggregate collections; and outbreak associations must be removed before a session can be deleted.



- Enter as much information as possible about the surveillance session you wish to delete. **Note:** Partial entry is acceptable for the **Session ID** and **Field Session ID**.

**Tip 1:** To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

**Tip 2:** To return to the previous screen without saving entered information, click **Cancel**.

**Tip 3:** To delete entered information, click **Clear**.

The screenshot shows the 'Vector Surveillance Session' search form. It includes fields for 'Session ID', 'Field Session ID', 'Status', 'Session Start Date Range' (From: 8/1/2019, To: 08/05/2019), 'Session Close Date Range' (From: , To: ), 'Vector Type', 'Species', 'Disease', 'Region' (Tbilisi), 'Rayon' (Didube-Chugureti), and 'Settlement'. A red circle with a '2' icon highlights the search criteria section. At the bottom, there are buttons for 'Cancel', 'Clear', 'Search', and 'Add'. A red circle with a '3' icon highlights the 'Search' button. A blue arrow with a 'T1' icon points to the 'Session Close Date Range' field.

- Click **Search**. The **Search Results** appear.

## Delete a Vector Surveillance Session

Pg. 2 of 2

**Tip 4:** To revise the **Search Criteria**, click the **Show/Hide** icon, make the necessary changes, and click **Search**.

**Tip 5:** To view additional information, click the **Show/Hide** icon next to the desired session.

- Click the **Session ID** corresponding to the desired session. The **Vector Surveillance Session** appears displaying all sections of the session.

**Tip 6:** To delete a detailed or aggregate collection, click the **Delete** icon corresponding to the desired collection, and click **Yes** on the dialog box that appears.

- Click **Delete**. A confirmation dialog box appears.

Vector Surveillance Session

**Search Criteria** T4

**Search Results**

Session ID	Status	Start Date	Close Date	Region	Rayon	Vector Type	
VSS190080	In Process	08/02/2019		Tbilisi	Didube-Chughureti	Tick	
VSS190081	In Process	08/01/2019		Tbilisi	Didube-Chughureti		
VSS190085	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190086	In Process	08/01/2019		Tbilisi	Didube-Chughureti		
VSS190090	In Process	08/01/2019		Tbilisi	Didube-Chughureti		

Number of records: 5 Page: 1

Cancel Add

Vector Surveillance Session

**Session Summary**

Session ID: VSS190072 Field Session ID: 23456t Status: In Process

**Location**

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Kakheti

Rayon: Telavi (rayon) Settlement: Zagarebi

Latitude (N.####): 0 Longitude (E.####): 0 Map:

Description of Location:

**Detailed Collections**

Add

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon	
<input type="checkbox"/>	1	VS190-032	Tick	7/10/2019	Kakheti	Telavi (rayon)	

Number of records: 1 Page: 1

**Aggregate Collections**

Add

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors	
2	Tick	7/4/2019	36	

Number of records: 1 Page: 1

Cancel Submit Delete

- Click **Yes**. A pop-up appears confirming successful deletion of the session.

Are you sure you want to delete this record?

6 Yes No