

Instructor Guide

Vector Module

EIDSS | Electronic Integrated Disease
Surveillance System

This page intentionally left blank



Module Overview

This module presents EIDSS functions related to entry and management of a vector surveillance session. To ensure pertinent information is shared across all sites, EIDSS records demographic information, geographical locations, laboratory analyses, sample tracking, epidemiological analyses, clinical information (including disease-specific clinical signs), and response measures.

This module consists of one lesson and its associated tasks:

- Vector Surveillance Session

Agenda

	Time
Module Overview Agenda	00:30
Lesson: Vector Surveillance Session Enter a Vector Surveillance Session and Detailed Collection Search for a Vector Surveillance Session Enter an Aggregate Collection Edit a Vector Surveillance Session Copy a Detailed Collection Delete a Vector Surveillance Session Exercise	01:30
Module Summary	00:30



Vector Surveillance Session

Overview

The user can enter, view, and edit the basic information about a vector surveillance session, including dates, location, vector types/species, samples, and laboratory tests. The data can be entered in detail, known as a detailed collection, or as a statistical report, known as an aggregate collection.

The following tasks are included in this lesson:

- Enter a Vector Surveillance Session and Detailed Collection
- Search for a Vector Surveillance Session
- Enter an Aggregate Collection
- Edit a Vector Surveillance Session
- Copy a Detailed Collection
- Delete a Vector Surveillance Session



Demonstration

Demonstrate each task for the students. After each demo, have them complete the task, inputting either the same data set used in the demonstration or different information. Allow them to refer to the quick guides as needed. Continue in this manner until all tasks are completed.



Exercise

Using the appropriate exercise have the students practice the tasks and apply the knowledge learned in this lesson. As students move through the exercise, they should record any IDs generated by the system for later use.

Note: There will be EIDSS fields for which information was not provided in the scenario. Tell students that they can leave these blank or add fictional, but realistic, data based on their epidemiological experience.



Enter a Vector Surveillance Session and Detailed Collection

Key Points


Include the task's key points in your demonstration.

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection Pg. 1 of 9

An EIDSS user with the appropriate permissions can create surveillance sessions in the database. This quick guide presents the steps to enter a vector surveillance session and detailed collection.

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears.



Tip 1: If there is concern that the desired session already exists in the database, perform a search before continuing.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

- Click **Add**. A blank **Vector Surveillance Session** screen appears.

3

- The **Field Session ID** is user generated, and the requirements may vary from country to country.
- The **Close Date** is only populated if the **Session Status** is set to **Closed**.
- The **Collection Effort** is the number of hours spent collecting vectors within the session.
- If **Exact Point** is selected as the **Location Type**, **Latitude** and **Longitude** will auto-populate based on the lowest address level chosen; the fields remain editable.
- If **Relative Point** is selected as the **Location Type** and all fields are completed, the system will auto-populate **Latitude** and **Longitude**; the fields are not editable.
- The **Distance** fields for **Relative Point** locations allow you to enter the number of kilometers or number of degrees from the chosen town or village.
- If the Vector Surveillance Session is created from an Outbreak Session, the outbreak location information auto-populates, but the fields are editable.

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 2 of 9

3. Complete the **Session Summary** and **Location** sections. **Note:** Fields marked with a red asterisk are mandatory. **Note:** Any commas used when entering **Longitude** and **Latitude** will appear as periods once the session is saved.

Tip 3: To navigate to the desired date, click the **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 4: To quickly enter **Latitude** and **Longitude**, click the **Map** icon and click the appropriate location on the map that appears (choosing a location on the map will overwrite any **Latitude** and **Longitude** already entered).

4. Click **Submit**. A window appears confirming successful entry of the session, and the screen refreshes to include the **Detailed Collections** and **Aggregate Collections** sections.
5. Click **OK**.

6

- The columns shown in the **Detailed Collection** grid can be customized by an administrator.

7

- Once a sample is added to the record, the **Vector Type** and **Species** are no longer editable.
- If **Exact Point**, **Relative Point**, or **Foreign Address** was selected as the **Type of Location** in the **Location** section, the same **Location Type** will be automatically selected for the **Collection Data** section. The information entered in the **Location** section will auto-populate; the fields remain editable.
- If **National** was selected in the **Location** section, **Exact Point** and **Relative Point** will be the only **Location Type** options in the **Collection Data** section.
- The **Host Reference** drop-down menu appears only if the chosen **Vector Type** is a pooled vector.
- Future dates are not allowed in any of the date fields.
- The **Collection Date** must be later than the **Start Date** and earlier than the **Identifying Date**, **Accession Date**, and **Result Date**.
- The **Identifying Date** must be later than the **Start Date** and **Collection Date** or earlier than the **Accession Date**, **Result Date**, and **Close Date**.

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 3 of 9

6. Click **Add** next to **Detailed Collections**. The **Vector Surveillance Session** screen appears displaying the **Collection Data** section.

Vector Surveillance Session

Session Summary

Session ID: VSS190085 Field Session ID: ASF-2019-43 Status: In Process

Detailed Collections

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon
No results were found.						

Number of records: Page:

Aggregate Collections

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors
No results were found.			

Number of records: Page:

Cancel Submit Delete

7. Complete the **Collection Data** section.

Tip 5: The **Session Summary** appears on all sections of the **Vector Surveillance Session** screen; to collapse and expand the **Session Summary**, click the **Show/Hide** icon.

8. Click **Next**. The **Vector Data** section appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190085 Field Session ID: ASF-2019-43 Status: In Process

Collection Data

Pool/Vector ID: Field Pool/Vector ID: ASF-2019-43-01

Vector Type: Flea

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Samegrelo-Zemo Svaneti

Rayon: Poti Settlement:

Latitude (#.####): Longitude (#.####): Map

Description of Location:

Elevation (m): 34 Surroundings: Forest - Deciduous trees

GEO Reference Source: Garmin eTrex Basis of Record: Preserved Specimen Host Reference:

Collected by Institution: NCDC Collector:

Collection Date: 8/10/2019 Collection Time: Day Collection Method: Adult, from open are Ectoparasites Collected:

Cancel Next

11

- This section is created from a flexible form and will vary significantly depending on the chosen vector.
- Choosing **Hare, Fox, Jackal, Wolf, or Other animals** as the **Vector Type** does not result in default data fields in the **Vector Specific Data** section.

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 4 of 9

9. Complete the **Vector Data** section.

Tip 6: The navigation menu allows you to access any section and displays a green checkmark next to each completed section.

Tip 7: To collapse the navigation menu, click the three lines at the top.

10. Click **Next**. The **Vector Specific Data** section appears.

11. Complete the **Vector Specific Data** section.

12. Click **Next**. The **Vector Surveillance Session** screen appears displaying all sections of the collection.

12

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 5 of 9

13. Review the information for accuracy.

Tip 8: To revise information, click the **Edit** icon in the appropriate section, make the necessary changes, and return to the **Review**.

14. Click **Submit**. A window appears confirming successful entry of the collection. **Note:** Once a detailed collection is saved the **Samples**, **Field Tests**, and **Laboratory Tests** sections appear.

Vector Surveillance Session

Session Summary

Session ID: VSS190085 Field Session ID: ASF-2019-43 Status: In Process

Collection Data

Pool/Vector ID: Field Pool/Vector ID: ASF-2019-43-01

Vector Type: Flea

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Samegrelo-Zemo Svaneti

Rayon: Poti Settlement:

Latitude (#.####): Longitude (#.####): Map

Description of Location:

Elevation (m): 34 Surroundings: Forest - Deciduous trees

GEO Reference Source: Basis of Record: Host Reference:

13

T8

Collection Data

Vector Data

Vector Specific Data

Review

Flea feeding status: Flea Flea life stage: Larva

Index of flea abundance:

Percent of occurrence:

Total stock of fleas on 1 hectare:

Cancel Submit

14

15. Click **OK**.

Vector Surveillance Session

☒ The record is saved successfully.

OK Return to Vector Surveillance Session Return to Dashboard

15

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 6 of 9

16. Click **Samples**. The **Samples** section appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190085 Field Session ID: ASF-2019-43 Status: In Process

Vector Data

Quantity: 32 Species: Amphipsylla georgica Sex: Mixed

Identified by Institution: Ministry of Agriculture Identified by:

Identifying Method: Morphology Identifying Date: 8/12/2019

Buttons: Cancel Previous Next Delete

Disclaimer | © 2016 - 2019 MoLHSA. All Rights Reserved. Build: 7.0.88.0 | Environment: Integration

Menu items: Collection Data, Vector Data, Vector Specific Data, **Samples**, Field Tests, Laboratory Tests, Review

17. Click **Add Sample**. The **Sample** window appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190085 Field Session ID: ASF-2019-43 Status: In Process

Samples

Add Sample

Bar	Code	Field Sample ID	Collection Date	Sample Type	Collected by	Sent To	Condition Received
No results were found.							

Buttons: Cancel Previous Next Delete

Disclaimer | © 2016 - 2019 MoLHSA. All Rights Reserved. Build: 7.0.88.0 | Environment: Integration

Menu items: Collection Data, Vector Data, Vector Specific Data, **Samples**, Field Tests, Laboratory Tests, Review

18

- The **Result Date** must be later than the **Collection Date**, **Identified Date**, and **Accession Date**.

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 7 of 9

18. Complete the **Sample** window.
19. Click **Save Sample**. **Note:** Repeat steps 17–19 to add additional samples.

Sample

Field Sample ID: A1-2019

Sample Type: Blood

Collection Date: 8/12/2019

Sent to Organization: Northrop Grumman

Collected by Institution: NCDC

Comment:

Accession Date:

Sample Condition Received:

Save Sample Close

20. Click **Next**. The **Field Tests** section appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190085

Field Session ID: ASF-2019-43

Status: In Process

Samples

Bar Code	Field Sample ID	Collection Date	Sample Type	Collected by	Sent To	Condition Received
SFB008190087	A1-2019	08/12/2019	Blood	NCDC	NGIT	

Number of records: 1

Page: 1

Cancel Previous Next Delete

20

Disclaimer | © 2016 - 2019 MoLHSA. All Rights Reserved. Build: 7.0.88.0 | Environment: Integration

21. Click **Add Field Test**. The **Field Test** window appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190085

Field Session ID: ASF-2019-43

Status: In Process

Field Tests

Add Field Test

Field Sample ID	Field Test	Test Category	Tested by Institution	Test Result	Test By	Diagnosis
No results were found.						

Cancel Previous Next Delete

21

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 8 of 9

22. Complete the **Field Test** window.
23. Click **Save Field Test**. **Note:** Repeat steps 21–23 to add additional field tests.

Field Test

Field Sample ID

Sample Type

Collection Date

Result Date

Test Name

Test Category

Tested by Institution

Tested By Person

Test Result

Disease

24. Click **Next**. The **Laboratory Tests** section appears. **Note:** This section is read-only.

Vector Surveillance Session

Session Summary

Session ID

Field Session ID

Status

Field Tests

Field ID	Sample	Field Test	Test Category	Tested By Institution	Test Result	Test By	Diagnosis
A1-2019	ASF Genotype	ASF Genotype			Gram positive		African Swine Fever

Number of records: 1 Page: 1

☒ Collection Data
☒ Vector Data
☒ Vector Specific Data
☒ Samples
☒ Field Tests
☒ Laboratory Tests
☐ Review

25. Click **Next**. The **Vector Surveillance Session** screen appears displaying all sections of the session.

Vector Surveillance Session

Session Summary

Session ID

Field Session ID

Status

Laboratory Tests

Lab Sample ID	Field Sample ID	Species	Test	Test Result	Disease	Sample Type	Field Collection Date	Test Category	Tested by Institution	Test By	Amendment History
No results were found.											

☒ Collection Data
☒ Vector Data
☒ Vector Specific Data
☒ Samples
☒ Field Tests
☒ Laboratory Tests
☐ Review

Instructor Notes

Enter a Vector Surveillance Session and Detailed Collection

Pg. 9 of 9

26. Review the information for accuracy.

27. Click **Submit**. A window appears confirming successful entry of the detailed collection.

Vector Surveillance Session

Session Summary

Session ID

VSS190085

Field Session ID

ASF-2019-43

Status

In Process

Collection Data

Pool/Vector ID

Field Pool/Vector ID

ASF-2019-43-01

Vector Type

Flea

Location Type

☒ Exact Point
 ☐ Relative Point
 ☐ Foreign Address
 ☐ National

Region

Samegrelo-Zemo Svaneti

Rayon

Poti

Settlement

Latitude (#.####o)

Longitude (#.####o)

Map

Description of Location

Elevation (m)

34

Surroundings

Forest - Deciduous trees

Host Reference

8/10/2019

Day

Adult, from open are

Samples

Bar Code	Field Sample ID	Collection Date	Sample Type	Collected by	Sent To	Condition Received
SFB008190087	A1-2019	08/12/2019	Blood	NCDC	NGIT	

Number of records: 1

Page: 1

Field Tests

Field ID	Sample	Field Test	Test Category	Tested By Institution	Test Result	Test By	Diagnosis
A1-2019	ASF Genotype				Gram positive		African Swine Fever

Number of records: 1

Page: 1

Laboratory Tests

Lab Sample ID	Field Sample ID	Species	Test	Test Result	Disease	Sample Type	Field Collection Date	Test Category	Tested by Institution	Test By	Amendment History
No results were found.											

Cancel

Submit

Delete



Search for a Vector Surveillance Session

Key Points

Include the task's key points in your demonstration.


Instructor Notes

Search for a Vector Surveillance Session

Pg. 1 of 1

An EIDSS user with the appropriate permissions can locate surveillance sessions in the database. This quick guide presents the steps to search for a vector surveillance session.

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears.



- Enter as much information as possible about the surveillance session for which you are searching. **Note:** Partial entry is acceptable for the **Session ID** and **Field Session ID**.

Tip 1: To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

Tip 3: To delete entered information, click **Clear**.

- Click **Search**. The **Search Results** appear displaying the vector surveillance sessions that match your **Search Criteria**.



Enter an Aggregate Collection

Key Points

Include the task's key points in your demonstration.


Instructor Notes

Enter an Aggregate Collection

Pg. 1 of 4

An EIDSS user with the appropriate permissions can add collections to surveillance sessions in the database. This quick guide presents the steps to enter an aggregate collection.

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears.



- Enter as much information as possible about the surveillance session to which you wish to add an aggregate collection. **Note:** Partial entry is acceptable for the **Session ID** and **Field Session ID**.

Tip 1: To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

Tip 3: To delete entered information, click **Clear**.

Vector Surveillance Session

Search Criteria

Use the following fields to search for existing session. Please fill in all known fields.

Session ID
VSS190084

Field Session ID

Status
▼

Session Start Date Range

From: [] To: [] **T1**

Session Close Date Range

From: [] To: []

Vector Type ▼ **Species** ▼ **Disease** ▼

Region ▼ **Rayon** ▼

Tbilisi Didube-Chugureti

Settlement ▼

Buttons: Cancel Clear Search Add

Annotations: T2 (Cancel), T3 (Clear), 3 (Search), 4 (Add)

- Click **Search**. The **Search Results** appear.

Instructor Notes

Enter an Aggregate Collection

Pg. 2 of 4

Tip 4: To revise the **Search Criteria**, click the **Show/Hide** icon, make the necessary changes, and click **Search**.

Tip 5: To view additional information, click the **Show/Hide** icon next to the desired session.

- Click the **Session ID** corresponding to the desired session. The **Vector Surveillance Session** screen appears displaying all sections of the session.
- Click **Add** next to **Aggregate Collections**. The **Vector Surveillance Session** screen appears displaying the **Aggregate Information** section.

Vector Surveillance Session

Search Criteria T4

Search Results

Session ID	Status	Start Date	Close Date	Region	Rayon	Vector Type
VSS190084	In Process	07/31/2019		Tbilisi	Didube-Chughureti	Tick

Number of records: 1 Page: 1

Cancel Add

Vector Surveillance Session

Session Summary

Session ID: VSS190084 Field Session ID: TST-8-12 Status: In Process

Location

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Tbilisi

Rayon: Didube-Chughureti Settlement:

Latitude (#####): 0 Longitude (#####): 0 Map

Description of Location:

Detailed Collections

Add

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon
<input type="checkbox"/>	17		Tick	8/12/2019	Tbilisi	Didube-Chughureti

Number of records: 1 Page: 1

Aggregate Collections

5 Add

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors
No results were found.			

Cancel Submit Delete

6

- The columns shown in the **Aggregate Collection** grid can be customized by an administrator.
- Changing the **Vector Type** will clear any data entered for the **Sex** and **Number of Pools/Vectors**.
- Once a **Disease** is added to the record, **Vector Type** and **Species** are no longer editable.
- The **Number of Positive Pools/Vectors** can be entered even if the **Number of Pools/Vectors** is left blank. If the **Number of Pools/Vectors** is entered, the **Number of Positive Pools/Vectors** for each disease cannot be greater than the number collected.

Instructor Notes

Enter an Aggregate Collection

Pg. 3 of 4

6. Complete the **Aggregate Information** section. **Note:** Fields marked with a red asterisk are mandatory.
7. Click **Submit**. A window appears confirming successful entry of the record and the screen refreshes displaying the **List of Diseases** section.

Vector Surveillance Session

Session Summary

Session ID: VSS190072 | Field Session ID: 23456t | Status: In Process

Aggregate Information

Record ID: [Text Field]

Collection Date: [Calendar Icon] | Vector Type: [Dropdown]

Species: [Dropdown] | Sex: [Dropdown] | Number of Pools/Vectors: [Dropdown]

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address

Region: [Dropdown] | Rayon: [Dropdown] | Settlement: [Dropdown]

Latitude (#.####): [Text Field] | Longitude (#.####): [Text Field] | Map: [Location Pin Icon]

Description of Location: [Text Field]

Cancel Submit

8. Click **OK**.

Vector Surveillance Session

The record is saved successfully. Record ID is VSS190084.

OK Return to Vector Surveillance Session Return to Dashboard

Instructor Notes

Enter an Aggregate Collection

Pg. 4 of 4

9. Choose a **Disease** from the drop-down menu.
10. Enter a **Number of Positive Pools/Vectors**.
11. Click **Add**. **Note:** Repeat steps 9–11 to add additional records.
12. Click **Submit**. A window appears confirming successful entry of the collection.

Vector Surveillance Session

Session Summary

Session ID

VSS190072

Field Session ID

23456t

Status

In Process

Aggregate Information

Record ID

VSS190072

Collection Date

11/5/2019

Vector Type

Tick

Species

Argasidae

Sex

Number of Pools/Vectors

Location Type

☒ Exact Point ☐ Relative Point ☐ Foreign Address

Region

Tbilisi

Rayon

Didube-Chugureti

Settlement

Latitude (#.####o)

Longitude (#.####o)

Map

Description of Location

List of Diseases

Disease

Lyme Disease

Number of Positive Pools/Vectors

4

Add

Number of records:

9

Page:

Cancel

Submit

Delete



Edit a Vector Surveillance Session

Explain that you will first demonstrate the steps to edit a session, and then demonstrate the steps to edit detailed and aggregate collections within that session.

Key Points

Include the task's key points in your demonstration.

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection


Pg. 1 of 6

An EIDSS user with the appropriate permissions can revise surveillance sessions in the database. This quick guide presents the steps to:

- Edit a vector surveillance session
- Edit a detailed collection
- Edit an aggregate collection

Edit a Vector Surveillance Session

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears.



- Enter as much information as possible about the surveillance session you wish to edit. **Note:** Partial entry is acceptable for the **Session ID** and **Field Session ID**.

Tip 1: To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

Tip 3: To delete entered information, click **Clear**.

The screenshot shows the 'Vector Surveillance Session' form. The 'Search Criteria' section is highlighted with a red box and a '2' label. It includes fields for 'Session ID', 'Field Session ID', 'Status' (dropdown), 'Session Start Date Range' (From/To), 'Session Close Date Range' (From/To), 'Vector Type', 'Species', 'Disease', 'Region', 'Rayon', and 'Settlement'. A 'T1' label points to the 'Session Start Date Range' field. At the bottom, there are buttons for 'Cancel', 'Clear', 'Search', and 'Add'. A 'T2' label points to the 'Cancel' button, a 'T3' label points to the 'Clear' button, and a '3' label points to the 'Search' button.

- Click **Search**. The **Search Results** appear.

5

- Changing the **Status** to **Closed** will cause the **Close Date** to auto-populate with the current date. Changing the **Status** back to **In Progress** will cause the **Close Date** to clear.
- The **Detailed Collections** and **Aggregate Collections** sections display any previously entered data and offer the ability to edit, delete, or add new detailed or aggregate collections.
- Changes made to the session **Location** will not affect the **Collection Data** location (if this location is different than the vector surveillance session **Location**).

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection

Pg. 2 of 6

Tip 4: To revise the **Search Criteria**, click the **Show/Hide** icon, make the necessary changes, and click **Search**.






Tip 5: To view additional information, click the **Show/Hide** icon next to the desired session.

- Click the **Session ID** corresponding to the desired session. The **Vector Surveillance Session** screen appears displaying all sections of the session.
- Make any necessary revisions to the **Session Summary** and **Location** sections. **Note:** The **Session Summary** defaults to collapsed; click the **Show/Hide** icon to expand the section.
- Click **Submit**. A window appears confirming successful modification of the session.

Vector Surveillance Session

Search Criteria T4

Search Results

Session ID	Status	Start Date	Close Date	Region	Rayon	Vector Type	
VSS190072	In Process	07/01/2019		Tbilisi	Didube-Chughureti	Flea	 T4
VSS190073	In Process	07/01/2019		Tbilisi	Didube-Chughureti	Fox	 T5
VSS190076	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190077	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190080	In Process	08/01/2019		Tbilisi	Didube-Chughureti		

Number of records: 6 Page: 1

Cancel Add

Vector Surveillance Session

Session Summary 5

Session ID: VSS190073 Field Session ID: Status: In Process

Outbreak ID:

Session Start Date: 7/1/2019

Description:


Collection Efforts:

Location

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region:

Rayon: Settlement:

Latitude (#.####): Longitude (#.####): Map: 

Description of Location:

Detailed Collections

Add

Cancel Submit 6

2

- Detailed collections, **Vector Type**, and **Species** cannot be edited in collections with associated samples.
- Sample data cannot be edited once a field test is associated.
- The **Sent to Organization** field cannot be changed if the associated sample has been accessioned (i.e., a **Lab Sample ID** was assigned).

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection

Pg. 3 of 6

Edit a Detailed Collection

1. Click the **Pool/Vector ID** corresponding to the detailed collection you wish to edit. The **Vector Surveillance Session** screen appears displaying all sections of the collection.

Vector Surveillance Session

Session Summary

Session ID: VSS190073 Field Session ID: Status: In Process

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon	
<input type="checkbox"/>	7		Fox	7/16/2019	Tbilisi	Didube-Chughureti	
<input type="checkbox"/>	8		Fox	7/15/2019	Tbilisi	Didube-Chughureti	

Number of records: 2 Page: 1

Aggregate Collections

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors
No results were found.			

Cancel Submit Delete

2. Click the **Edit** icon in the appropriate section. The chosen section appears.

Vector Surveillance Session

Session Summary

Session ID: VSS190073 Field Session ID: Status: In Process

Collection Data

Pool/Vector ID: Field Pool/Vector ID: Vector Type: Fox

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Tbilisi Rayon: Didube-Chughureti Settlement: Latitude: Longitude: Map

Description of Location:

Laboratory Tests

Lab Sample ID	Field Sample ID	Species	Test	Test Result	Disease	Sample Type	Field Collection Date	Test Category	Tested by Institution	Test By	Amendment History
---------------	-----------------	---------	------	-------------	---------	-------------	-----------------------	---------------	-----------------------	---------	-------------------

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection

Pg. 4 of 6

3. Make any necessary revisions.
4. Use the navigation menu to access and revise other sections as needed.
5. Click **Review**. The **Vector Surveillance Session** screen appears displaying all sections of the detailed collection.

Vector Surveillance Session

Session Summary

Session ID: VSS190073 Field Session ID: Status: In Process

Collection Data

Pool/Vector ID: Field Pool/Vector ID: **3**

Vector Type: Fox

Location Type: ☒ Exact ☐ Relative ☐ Foreign ☐ National

Region: Tbilisi

Rayon: Didube-Chugureti Settlement: **4**

Latitude (#.####): Longitude (#.####): Map

Description of Location:

Elevation (m): 0 Surroundings: **5**

GEO Reference Source: Basis of Record: Host Reference:

Collected by Institution: Kakheti Region Veterinary Office Collector: **5**

Collection Date: 7/16/2019 Collection Time: Collection Method: Ectoparasites Collected: **5**

Navigation Menu:

- Collection Data **4**
- Vector Data
- Vector Specific Data
- Samples
- Field Tests
- Laboratory Tests
- Review**

Cancel Next Delete

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection

Pg. 5 of 6

6. Review the revised information for accuracy.
7. Click **Submit**. A window appears confirming successful modification of the detailed collection.

Vector Surveillance Session

Session Summary

Session ID

VSS190073

Field Session ID

Status

In Process

Collection Data

Pool/Vector ID

Field Pool/Vector ID

Vector Type

Fox

Location Type

☒ Exact Point
 ☐ Relative Point
 ☐ Foreign Address
 ☐ National

Region

Tbilisi

Rayon

Didube-Chugureti

Settlement

Latitude (#.####o)

Longitude (#.####o)

Map

Description of Location

Elevation (m)

0

Surroundings

GEO Reference Source

Basis of Record

Host Reference

Collected by Institution

Kakheti Region Veterinary Office

Collector

Collection Date

7/16/2019

Collection Time

Day

Collection Method

By hand

Ectoparasites Collected

Vector Data

Quantity

0

Species

Vulpes vulpes

Sex

Identified by Institution

Identified by

Field Sample ID

Field Test

Test Category

Tested By Office

Test Result

Test By

Diagnosis

No results were found.

Laboratory Tests

Lab Sample ID

Field Sample ID

Species

Test

Test Result

Disease

Sample Type

Field Collection Date

Test Category

Tested by Institution

Test By

Amendment History

No results were found.

Cancel

Submit

Delete

7

Instructor Notes

Edit a Vector Surveillance Session, Detailed Collection, and Aggregate Collection

Pg. 6 of 6

Edit an Aggregate Collection

1. Click the **Record ID** corresponding to the aggregate collection you wish to edit. The **Vector Surveillance Session** screen appears displaying the collection.

Vector Surveillance Session

Session Summary

Session ID: VSS190073 Field Session ID: Status: In Process

Location

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Tbilisi

Rayon: Didube-Chughureti Settlement:

Latitude (N.####c): 0 Longitude (E.####c): 0 Map

Description of Location:

Detailed Collections

Copy Pool/Vector ID Field Pool/Vector ID Vector Type Collection Date Region Rayon Add

<input type="checkbox"/>	7		Fox	7/16/2019	Tbilisi	Didube-Chughureti			
<input type="checkbox"/>	8		Fox	7/15/2019	Tbilisi	Didube-Chughureti			

Number of records: 2 Page: 1

Aggregate Collections

Record ID Vector Type Collection Date Number of Positive Pools/Vectors Add

14	Flea	8/7/2019	45			
----	------	----------	----	--	--	--

Number of records: 1 Page: 1

2. Make any necessary revisions.
3. Click **Submit**. A window appears confirming successful modification of the aggregate collection.

Vector Surveillance Session

Session Summary

Session ID: VSS190073 Field Session ID: Status: In Process

Aggregate Information

Record ID: VSS190073

Collection Date: 8/7/2019 Vector Type: Flea

List of Diseases

Disease: Number of Positive Pools/Vectors: Add

Diagnosis	Quantity
African Swine Fever	3

Number of records: 1 Page: 1

Cancel Submit



Copy a Detailed Collection

The ability to copy is available for detailed collections only.

Key Points

Include the task's key points in your demonstration.


Instructor Notes

Copy a Detailed Collection

Pg. 1 of 3

An EIDSS user with the appropriate permissions can create duplicates of detailed collections and modify them as needed in the database. This quick guide presents the steps to copy a detailed collection.

1. From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears.



2. Enter as much information as possible about the surveillance session with the detailed collection you wish to copy.
Note: Partial entry is acceptable for the **Session ID** and **Field Session ID**.

Tip 1: To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

Tip 3: To delete entered information, click **Clear**.

3. Click **Search**. The **Search Results** appear.

The screenshot shows the 'Vector Surveillance Session' form. It has a 'Search Criteria' section with the following fields: Session ID, Field Session ID, Status (dropdown), Session Start Date Range (From and To date pickers), Session Close Date Range (From and To date pickers), Vector Type (dropdown), Species (dropdown), Disease (dropdown), Region (dropdown), Rayon (dropdown), and Settlement (dropdown). A red circle with the number '2' points to the Session ID field. At the bottom of the form, there are four buttons: Cancel, Clear, Search, and Add. A red circle with the number '3' points to the Search button. There are also blue callout boxes: 'T1' points to the Session Close Date Range To field, 'T2' points to the Cancel button, and 'T3' points to the Clear button.

Instructor Notes

Copy a Detailed Collection

Pg. 2 of 3

Tip 4: To revise the **Search Criteria**, click the **Show/Hide** icon, make the necessary changes, and click **Search**.






Tip 5: To view additional information, click the **Show/Hide** icon next to the desired session.

- Click the **Session ID** corresponding to the desired session. The **Vector Surveillance Session** screen appears displaying all sections of the session.
- Check the **Copy** box of the desired detailed collection. The **Copy Detailed Collection Record** window appears.

Vector Surveillance Session

Search Criteria T4

Search Results

Session ID	Status	Start Date	Close Date	Region	Rayon	Vector Type	
VSS190072	In Process	07/01/2019		Tbilisi	Didube-Chughureti	Flea, Hare	
VSS190073	In Process	07/01/2019		Tbilisi	Didube-Chughureti	Flea, Fox	
VSS190076	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190077	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190080	In Process	08/01/2019		Tbilisi	Didube-Chughureti	Tick	

Number of records: 6 Page: 1

Cancel Add

Vector Surveillance Session

Session Summary

Session ID: VSS190072 Field Session ID: Status: In Process

Location


Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Tbilisi Rayon: Didube-Chughureti Settlement:

Latitude (#.####): Longitude (#.####): Map



Description of Location:

Detailed Collections

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon	
<input type="checkbox"/>	9		Flea	7/22/2019	Tbilisi	Didube-Chughureti	

Number of records: 1 Page: 1

Aggregate Collections

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors	
13	Flea	7/22/2019	5	
15	Hare	8/7/2019	90	

Number of records: 2 Page: 1

Cancel Submit Delete

6

- The **Collection Data** box is automatically checked and cannot be unchecked.
- Checking the **Vector Data** box will copy all data from the **Vector Data** section of the parent record to the child record.
- Checking the **Sample** box will copy data in the **Sample Type, Collection Date, Collected By Institution, and Sent to Organization** fields of the parent record to the child record. All copied data is editable.
- Checking the **Field Tests** box will copy data in the **Sample Type, Collection Date, Test Name, Test Category, Tested by Institution, and Tested by Person** fields from the parent record to the child record. All copied data is editable. If the **Field Tests** section is copied, the **Samples** section will also be copied.

Instructor Notes

Copy a Detailed Collection

Pg. 3 of 3

6. Check the box of each section you wish to copy. **Note:** The **Collection Data** section is copied automatically.
7. Click **Submit**. The **Vector Surveillance Session** screen appears displaying the copied record in the **Detailed Collections** section with a new **Pool/Vector ID** assigned.

Copy Detailed Collection Record

Detailed Collection Id:6

☒ Collection Data

☒ Vector Data

☒ Samples

☒ Field Tests

Submit Cancel



Delete a Vector Surveillance Session

A vector surveillance session cannot be deleted until any associated aggregate or detailed collections are deleted.

Key Points

Include the task's key points in your demonstration.


Instructor Notes

Delete a Vector Surveillance Session

Pg. 1 of 2

An EIDSS user with the appropriate permissions can remove surveillance sessions from the database. This quick guide presents the steps to delete a vector surveillance session.

- From the **EIDSS Dashboard**, choose **SURVEILLANCE SESSION** from the **VECTOR** drop-down menu. The **Vector Surveillance Session** screen appears. **Note:** All field tests and samples within detailed collections; detailed collections; aggregate collections; and outbreak associations must be removed before a session can be deleted.



- Enter as much information as possible about the surveillance session you wish to delete. **Note:** Partial entry is acceptable for the **Session ID** and **Field Session ID**.

Tip 1: To navigate to the desired date, click a **Calendar** icon and then the **Select Month > Select Year > Select Decade** bars at the top of the calendar.

Tip 2: To return to the previous screen without saving entered information, click **Cancel**.

Tip 3: To delete entered information, click **Clear**.

The screenshot shows the 'Vector Surveillance Session' search criteria form. The form is titled 'Search Criteria' and includes a note: 'Use the following fields to search for existing session. Please fill in all known fields.' The form contains several input fields and dropdown menus, including 'Session ID', 'Field Session ID', 'Status', 'Session Start Date Range' (From: 8/1/2019, To: 08/05/2019), 'Session Close Date Range' (From: , To:), 'Vector Type', 'Species', 'Disease', 'Region' (Tbilisi), 'Rayon' (Didube-Chugureti), and 'Settlement'. A red circle with a '2' icon highlights the 'Field Session ID' field. A blue arrow with a 'T1' icon points to the 'Session Close Date Range' 'To' field. At the bottom, there are buttons for 'Cancel', 'Clear', 'Search', and 'Add'. A red circle with a '3' icon highlights the 'Search' button. A blue arrow with a 'T2' icon points to the 'Cancel' button, and a blue arrow with a 'T3' icon points to the 'Clear' button.

- Click **Search**. The **Search Results** appear.

5

- Deleting a sample from a detailed collection will remove it from the surveillance session, but it will remain in the Laboratory module.

Instructor Notes

Delete a Vector Surveillance Session

Pg. 2 of 2

Tip 4: To revise the **Search Criteria**, click the **Show/Hide** icon, make the necessary changes, and click **Search**.

Tip 5: To view additional information, click the **Show/Hide** icon next to the desired session.

- Click the **Session ID** corresponding to the desired session. The **Vector Surveillance Session** appears displaying all sections of the session.

Tip 6: To delete a detailed or aggregate collection, click the **Delete** icon corresponding to the desired collection, and click **Yes** on the dialog box that appears.

- Click **Delete**. A confirmation dialog box appears.

Vector Surveillance Session

Search Criteria T4

Search Results

Session ID	Status	Start Date	Close Date	Region	Rayon	Vector Type	
VSS190080	In Process	08/02/2019		Tbilisi	Didube-Chughureti	Tick	
VSS190081	In Process	08/01/2019		Tbilisi	Didube-Chughureti		
VSS190085	In Process	08/05/2019		Tbilisi	Didube-Chughureti		
VSS190086	In Process	08/01/2019		Tbilisi	Didube-Chughureti		
VSS190090	In Process	08/01/2019		Tbilisi	Didube-Chughureti		

Number of records: 5 Page: 1

Cancel Add

Vector Surveillance Session

Session Summary

Session ID: VSS190072 Field Session ID: 23456t Status: In Process

Location

Location Type: ☒ Exact Point ☐ Relative Point ☐ Foreign Address ☐ National

Region: Kakheti

Rayon: Telavi (rayon) Settlement: Zagarebi

Latitude (N.####): 0 Longitude (E.####): 0 Map:

Description of Location:

Detailed Collections

Add

Copy	Pool/Vector ID	Field Pool/Vector ID	Vector Type	Collection Date	Region	Rayon	
<input type="checkbox"/>	1	VS190-032	Tick	7/10/2019	Kakheti	Telavi (rayon)	

Number of records: 1 Page: 1

Aggregate Collections

Add

Record ID	Vector Type	Collection Date	Number of Positive Pools/Vectors	
2	Tick	7/4/2019	36	

Number of records: 1 Page: 1

Cancel Submit Delete

- Click **Yes**. A pop-up appears confirming successful deletion of the session.

Are you sure you want to delete this record?

6 Yes No



Module Summary

This module focused on EIDSS functions related to entering and managing a vector surveillance session. The lesson was presented along with its associated tasks. Students practiced these tasks after instructor demonstrations and then applied the knowledge they learned in an exercise at the end of the lesson.